

INTRODUCING CAPP Expell & ExpellPlus **TIPS and TUBES**

Capp TIPS are manufactured in Germany in clean rooms and certified to be RNase, DNase, DNA and Pyrogen free. Capp offers a wide portfolio of *Bulk, Reload and Filter Low retention* tips, universally compatible with most pipette brands.

Buy any 3
boxes of
tips and tubes

Get a 4th
box of your
choice

FREE



Bulk Tips

5030010	10 µL Sterile Tips Clear	×1,000
5130010	10 µL Clear Tips	1,000
5130040	10 µL XL Clear Tips	1,000
5130070	200 µL Clear Tips	1,000
5130100	300 µL Clear tips	1,000
5130130	1,000 µL Clear Tips	1,000
4130075	200 µL Yellow Tips	1,000
4130135	1,000 µL Blue Tips	1,000
5130160	5,000 µL Clear Tips	250
5130190	10,000 µL Clear Tips	200

Filter Low Retention Tips

5030030	10 µL Sterile Filter Tips	10×96
5030060	10 µL XL Sterile Filter Tips	10×96
5030062	20 µL Sterile Filter Tips	10×96
5030066	100 µL Sterile Filter Tips	10×96
5030090	200 µL Sterile Filter Tips	10×96
5030115	300 µL Sterile Filter Tips	10×96
5130150	1,000 µL Sterile Filter Tips	8×96
5130123	1,250 µL Sterile Filter Tips	8×96

Reload Tips

50030035	10 µL Reload Tips	10×96
5030095	200 µL Reload Tips	10×96
5130155	1,000 µL (1,250 µL) Reload Tips	5×96

For more information,
prices or samples,
please visit:
pocdscientific.com.au

The world's best medical technology — in the lab or at the point of care

Order online at www.pocdscientific.com.au or call 1800 640 075

POCD Scientific | Unit B, 19-21 Loyalty Road, North Rocks NSW 2151 | info@pocdscientific.com.au

CAPP Expell Microcentrifuge Tubes

Manufactured in Germany, Capp Tubes are made of virgin Polypropylene and offer perfect sealing and evaporation protection.

All Expell Microcentrifuge tubes are DNase, RNase, Pyrogen free and can withstand up to 20,000 g.

Polished interior surface prevents Protein and Nucleic acid adhesion.

High clarity for good sample visibility.




Range to choose from:

PART No.	Description	Pack Size
5100500	Microcentrifuge tubes 0.5 mL	1,000 pk
5001500	Microcentrifuge tubes Low Binding 1.5 mL	500 pk
5101500	Microcentrifuge tubes Secure Safelock 1.5 mL	500 pk
5102000	Microcentrifuge tubes 2 mL	500 pk

Tips to choose from:

- [Capp Bulk Standard](#)
- [Capp Low Retention](#)
- [Capp Reloads](#)



For more information, a demo, or to place an order, please call:

NSW 0418 231 583 | VIC 0418 270 636 | QLD 0413 590 215

The world's best medical technology — in the lab or at the point of care

Order online at www.pocdscientific.com.au or call **1800 640 075**

POCD Scientific | Unit B, 19-21 Loyalty Road, North Rocks NSW 2151 | info@pocdscientific.com.au